

## Western Blot

1. Thaw sample on the ice.
2. Turn on the heater at 95°C.
3. Make gel with cassette (in vitro) and prepare samples.
  - 1) Separating gel: 7.5% polymerase gel with 4% shocker. Receipt:
    - dd-water 4.14 ml
    - Tris HCL 1.875 ml PH=8.8
    - 20% SDS 37.5 ul
    - A/Bis 1.405 ml
    - 10% APS 37.5 ul
    - TEMED 5 ulPipette the gel into cassette until it reaches 2<sup>nd</sup> line.  
Add water sol. Bustinol (not too much) to prevent evaporation.  
Wait 30 min
  - 2) Receipt: tissue, lysis buffer (stored at -20°C), loading dye (stored at -20°C)  
Volume: 25 ul, Protein: 50 ug  
Then, add 1/20 volume  $\beta$ -ME  
Centrifuge for 1 min, heat 5 min or more at 95°C, stored at RT.
  - 3) Discard the Water sol. Bastinol, rinse thoroughly with dd-water, wipe off.  
Loading gel: 4% polymerase gel. Receipt:
    - dd-water 1.92 ml (RT)
    - Tris HCL 0.75 ml (RT) PH=6.8
    - 20% SDS 15 ul (RT)
    - A/Bis 0.3 ml (4°C)
    - 10% APS 15 ul (-20°C)
    - TEMED 2 ul (RT)
  - 4) Move the comb and rinse thoroughly.
4. Load 1\* SDS Elution buffer between the gels until it reaches the first line.  
Receipt: 10 \* SDS elution buffer stock
  - 30.3 g Tris base (0.25M)
  - 144 g Glycine (1.92M)
  - 1.0 g SDSDilute to 1000 ml with dd-water  
Stored at 4 °C  
1:10 diluted in dd-water when used.
5. Load the marker (ladder) and samples.
6. Load 1\* SDS Elution buffer outside the gels until the samples were immersed
7. Run the gel
  - 75V 5~10 min until blue dye pass the first line
  - 120 V 1.5 hr -> 2.5 hr until blue dye is approximately 0.5-1 mm above the bottom.Note: if the blue dye is not straight, the gel does not perfectly polymerase
8. Transfer the protein to membrane.  
Pre-wet the sponge and filter paper, membrane with 1/2 \* Transfer buffer stock.
  - Sponge \*2 (blotting pad)
  - Filter paper
  - Transfer membrane

- Gel
- Filter paper
- Sponge \*2 (blotting pad)

Put a sponge between them if transferring two membrane at the same time.

And run overnight at 15 V in cold room.

Receipt: 10 X transfer buffer

30.3 g Tris Base (0.25 M)

144 g Glycine (1.92 M)

Diluted to 1000 ml with dd-water PH=8.2-8.3. stored at 4°C

½ \* Transfer buffer

50 ml 10\* transfer buffer

50 ml Methanon (RT Hood)

900 ml water Stored at 4°C

9. remove the membrane, rinse with TBST (Tris -buffered Tween-20)

Receipt: TBST

100 ml TBS 10X

900 ml dd-water

500 ul tween-20 ( hood)

10X TBS (Tris buffered saline)

24.2 g Trizma (0.2m)

80 g sodium chloride (1.37M)

PH to 7.6 using 6M HCL

Diluted to 1000 ml with dd-water stored at RT

10. Wash membrane in TBST

11. wash membrane in TBSTM 1 hour to block.

Recipe: Tris buffered saline tween -20 w milk

10 g dry milk

20 ml 10X TBS

100 ul Tween- 20

180 ml dd H<sub>2</sub>O

12. Incubated in Primary antibody 1 hour

13. wash membrane in TBST 4 times. 5 min/ each

14. Incubated in secondary antibody 2 hours

15. Wash membrane with TBST 4 times. 5min/each

16. Turn on film developer (standby), prepare ECL reagent. (5ml, 5ml)

17. Expose for 2 min – 5 min, and run